

Editing Tools

The screenshot shows the Flare software interface with several editing tools and shortcuts highlighted:

- Top Bar:** Ctrl+S, Ctrl+Z, Ctrl+Y, Ctrl+V, Ctrl+C, Ctrl+X, F5.
- Editing Tools:** Delete, Edit Residue, Fix Residue Order, Del H, Charge for pH7, Add H, Add Fields, Optimize Alignment, Lasso, Paste, Copy Picked, Cut Picked, Clipboard.
- Ligands Panel:** Lines (L), Thin Stick (T), Capped Stick (C), Ball and Stick (B), CPK (K).
- Annotations:**
 - Del, Backspace (delete all picked atoms)
 - Hold Ctrl to toggle Lasso mode temporarily
 - Hover mouse to change bond order: type 1/2/3
 - LMB-click to change bond order: 1/2/3
 - LMB-drag to rotate smaller end, Shift+LMB-drag to rotate bigger end
 - LMB-drag to sprout a new atom or ring
 - Hover mouse to change charge status: type +, -,
 - Hover mouse to change element: type N/O/S/F..

Shortcut	Action
</> keys	Zoom in/out
RMB-click on any atom	Display long right-click menu
Ctrl+LMB	Lasso to select part or all of the ligand
Home key	Center view
+,/,-, _	Add positive/negative ligand field surfaces

Shortcut	Action
LMB-drag on white space	Rotate the 3D view
Shift+LMB-drag on white space	Rotate molecule/selection relative to other objects
RMB-drag on white space	Translate the 3D view
Shift+RMB-drag on white space	Translate molecule relative to other objects
Shift+RMB-drag on any atom of selection	Translate a selection of the molecule

Shortcut	Action
Pick water molecule, then Shift+LMB-drag on oxygen	Reorient water molecule
LMB-click on atom + atom selected in drawing widget	Change atom type
LMB-click on atom + ring selected in drawing widget	Grow a ring
RMB-drag on atom	Move selected atom
Ctrl+A	Select all molecule